

Fine quality, Integrity and service

Aluminum wheelchair series



1. Flat Tube Design, increase compression capacity



2. PU Rear Wheel with 36 Steel Wire, more durable



3. Manual Lock, humanized and safer

Fixed PU armrest and footrest. Aluminum drop back handle. Reticulated and breathable cushion. 6" PU front wheel and 22" PU rear wheel. With hand brake.



KJT107A >>



Model No	67	34	44	56	15	48	91	100	40	46	100 (Kg)	DIM 95x33x73	NW GW 12.2 13.7	CNTS 20' GP FCL 148	CNTS 40' HQ FCL 289
KJT107A	67	34	44	56	15	48	91	100	40	46	100	95x33x73	12.2 13.7	148	289



KJT107B >>



Thicken the Tube Wall



Aluminum Drop Back Handle

Fixed PU armrest and footrest. Aluminum drop back handle. Reticulated and breathable cushion. 6" PU front wheel and 12" PU rear wheel. With hand brake.



KJT107C >>



Brake Dust Cover



Fashion Foot-tape Design

Fixed PU armrest and footrest. Reticulated and breathable cushion. 6" front wheel and 22" rear wheel with solid tire. With hand brake.

Model No	66	31	44	30	15	46	87	97	40	46	100 (Kg)	DIM 84x31x73	NW GW 11.5 13.5	CNTS 20' GP FCL 143	CNTS 40' HQ FCL 347
KJT107B	66	31	44	30	15	46	87	97	40	46	100	84x31x73	11.5 13.5	143	347
KJT107C	67	31	44	56	15	46	85	94	40	46	100	96x33x88	12.5 14.5	97	237



★ FY863LF1

Aluminum Frame, Drop Back Handle with Locking United Brake, Fixed Inclined Armrest, Flip Up Footrest, Solid Castor, Rear MAG Wheel with PU Tyre, Plastic Side Board, Plastic Footplate.
铝合金车架, 后折把手连刹车(带自锁功能), 固定斜扶手, 上翻脚托, 实心前轮, 塑料铝框PU后轮, 塑料护板, 塑料脚踏板



★ FY867LF1

Aluminum Frame, Drop Back Handle with Locking United Brake, Fixed Inclined Armrest, Flip Up Footrest, Solid Castor, Rear Spoke Wheel with PU Tyre, Plastic Hand Rim, Plastic Side Board, Plastic Footplate.
铝合金车架, 后折把手连刹车(带自锁功能), 固定斜扶手, 上翻脚托, 实心前轮, 钢丝铝框PU后轮, 塑料手轮, 塑料护板, 塑料脚踏板



★ FY872L

Aluminum Frame with Locking United Brake, Fixed Armrest & Footrest, Solid Castor, Rear MAG Wheel with PU Tyre, Plastic Footplate.
铝合金车架, 后把连刹车(带自锁功能), 固定扶手, 固定脚托, 实心前轮, 塑料铝框PU后轮, 塑料脚踏板



★ FY871L

Aluminum Frame with Locking United Brake, Fixed Armrest & Footrest, Solid Castor, Rear Spoke Wheel with PU Tyre, Plastic Hand Rim & Footplate.
铝合金车架, 后把连刹车(带自锁功能), 固定扶手, 固定脚托, 实心前轮, 钢丝铝框PU后轮, 塑料手轮, 塑料脚踏板



★ FY977LF1

Aluminum Frame, Drop Back Handle with Locking United Brake, Flip Up & Detachable Armrest, Detachable Footrest, PU Castor, Rear MAG Wheel with PU Tyre, Plastic Side Board & Footplate.
铝合金车架, 后折把手连刹车(带自锁功能), 上翻可拆扶手, 可拆脚托, PU前轮, 塑料铝框PU后轮, 塑料护板, 塑料脚踏板



★ FY977LF4

Aluminum Frame, Flip Up Inclined Armrest, Detachable Footrest, Solid Castor, Rear Spoke Wheel with PU Tyre, Plastic Hand Rim & Side Board, Alum. Footplate.
铝合金车架, 上翻斜扶手, 可拆脚托, 实心前轮, 钢丝铝框PU后轮, 塑料手轮, 塑料护板, 塑料脚踏板



★ FY864LBP

Aluminum Frame, Fixed Armrest, Fixed Footrest, Solid Castor, Rear MAG Wheel with PU Tyre, Plastic Hand Rim, Side Board & Footplate.
铝合金车架, 固定扶手, 固定脚托, 实心前轮, 塑料铝框PU后轮, 塑料手轮, 塑料护板, 塑料脚踏板



★ FY864L

Aluminum Frame, Fixed Armrest, Fixed Footrest, Solid Castor, Rear Spoke Wheel with PU Tyre, Plastic Side Board, Aluminum Footplate.
铝合金车架, 固定扶手, 固定脚托, 实心前轮, 钢丝铝框PU后轮, 塑料护板, 铝合金脚踏板

MAIN PARAMETERS (CM) 主要参数(单位:厘米)

MODEL NO.	座宽	座深	座高	座宽	座高	脚踏高	脚踏宽	脚踏深	脚踏高	脚踏宽	脚踏深	脚踏高	脚踏宽	脚踏深
FY863LF1	57	46	38	15	49	88	100	74*27*72	10.00	11.50	164	452		
FY867LF1	82	46	51	15	49	88	100	79*27*72	11.00	12.50	175	420		
FY977LF1	54	41	36	15	50	91	100	70*26*73	12.80	14.80	192	485		
FY977LF4	61	44	51	15	49.5	91	100	79*25*91	12.20	13.70	154	360		

MAIN PARAMETERS (CM) 主要参数(单位:厘米)

MODEL NO.	座宽	座深	座高	座宽	座高	脚踏高	脚踏宽	脚踏深	脚踏高	脚踏宽	脚踏深	脚踏高	脚踏宽	脚踏深
FY872L	63	46	30	15	47	86	100	86*26*68	09.00	11.50	124	318		
FY871L	63	46	51	15	47	86	100	86*27*68	10.50	13.00	124	318		
FY864LBP	66	46	60	20	48	88	100	99*28*69	11.80	14.00	108	272		
FY864L	65	46	60	20	48	88	100	99*28*69	11.80	14.00	108	272		



★ FY976LF1

Aluminum Frame with Locking United Brake, Flip Up Armrest & Footrest, Solid Castor, Rear Solid Wheel.

铝合金车架, 后把连刹车(带自锁功能), 上翻扶手, 上翻脚踏, 实心前轮, 实心后轮, 轻便可折叠, 免托运上飞机



★ FY976LF2

Aluminum Frame with Locking United Brake, Flip Up Armrest, Plastic Footplate, Solid Castor, Rear Solid Wheel.

铝合金车架, 后把连刹车(带自锁功能), 上翻扶手, 塑料脚踏板, 实心前轮, 实心后轮, 轻便可折叠, 免托运上飞机



★ FY802

Powder Coating Steel Frame, Drop Back Handle with Locking United Brake, Fixed Inclined Armrest, Flip Up Footrest, Solid Castor, Rear MAG Wheel with PU Tyre, Plastic Side Board & Footplate.

喷涂钢制车架, 后折把手连刹车(带自锁功能), 固定斜扶手, 上翻脚踏, 实心前轮, 塑料轮毂PU后轮, 塑料护板, 塑料脚踏板



★ FY802F1

Powder Coating Steel Frame, Drop Back Handle with Locking United Brake, Fixed Inclined Armrest, Flip Up Footrest, Solid Castor, Rear Spoke Wheel with PU Tyre, Plastic Side Board, Hand Rim & Footplate.

喷涂钢制车架, 后折把手连刹车(带自锁功能), 固定斜扶手, 上翻脚踏, 实心前轮, 钢丝铝框PU后轮, 塑料护板, 塑料手轮, 塑料脚踏板



★ FY869L

Aluminum Frame, Fixed Inclined Armrest, Fixed Footrest, Solid Castor, Rear Spoke Wheel with PU Tyre, Plastic Hand Rim, Side Board & Footplate

铝合金车架, 固定斜扶手, 固定脚踏, 实心前轮, 钢丝铝框PU后轮, 塑料手轮, 塑料护板, 塑料脚踏板(标配: 格布)



★ FY871LABJ

Aluminum Frame, Drop Back Handle with Locking United Brake, Fixed Armrest & Footrest, Solid Castor, Rear MAG Wheel with PU Tyre, Plastic Hand Rim & Footplate.

铝合金车架, 后折把手连刹车(带自锁功能), 固定扶手, 固定脚踏, 实心前轮, 塑料轮毂PU后轮, 塑料手轮, 塑料脚踏板



★ FY809E

1x40"HQ CTNR LOADING 570PCS. Chromed Steel Frame, Detachable Pushing Handle, Fixed Armrest, Detachable Footrest, Solid Castor, Solid Rear Wheel, Aluminum Footplate.

电镀钢制车架, 可拆推手把, 固定扶手, 可拆脚踏, 实心前轮, 实心后轮, 铝合金脚踏板, 高柜可装570台



★ FY809EF1

Powder Coating Steel Frame, Fixed Armrest, Detachable Footrest, Solid Castor, Solid Rear Wheel, Plastic Hand Rim & Footplate.

喷涂钢制车架, 固定扶手, 可拆脚踏, 实心前轮, 实心后轮, 塑料手轮, 塑料脚踏板

MAIN PARAMETERS (CM) 主要参数(单位:厘米)

MODEL NO.	座宽	座深	座高	座宽	座高	座深	座高	座宽	座高	座深	座高	座宽	座高	座深	座高	座宽	座高	座深	座高	
FY976LF1	51	43	19	15	50	94	75	52*25*73	07.00	08.00	320	720								
FY976LF2	51	40	15	10	48	91	75	52*28*61	07.00	08.40	336	776								
FY859L	82	44	60	20	50	88	100	100*28*68	11.20	13.20	84	248								
FY871LABJ	63	46	51	15	47	86	100	86*27*88	10.50	13.00	124	316								

MAIN PARAMETERS (CM) 主要参数(单位:厘米)

MODEL NO.	座宽	座深	座高	座宽	座高	座深	座高	座宽	座高	座深	座高	座宽	座高	座深	座高	座宽	座高	座深	座高	
FY802	62	46	30	15	47	81	75	71*31*71	11.40	13.40	166	425								
FY802F1	64	46	51	15	47	81	75	76*33*71	12.00	14.00	150	378								
FY809E	63	46	58	20	47	86	100	77*21*81	16.50	17.50	200	570								
FY809EF1	63	46	58	20	49	87	100	78*21*88	16.50	17.50	200	570								