

KN95 MEDICAL FACE MASK



Hidden soft nose clip: can follow facial contour adjustment, fit the face.

Elastic earloop with dural point attachment helps provide a secure seal.

Disposable, single-use, breathable, light and environment-friendly.

PRODUCT PARAMETERS

Specification

Nose Clip: Made of plastic materials/Iron Wire
Material: Melt-blown Fabric+Hot Air Cotton + Non-woven Fabric
Sterile: None
Flammability: No Afterflame
Supply Ability: 500000 Piece/Pieces per Day
Standard: GB2626-2006 Certificates: ISO

Type: Earloop
Shape: Flat-Fold
BEF: $\geq 95\%$
Filtration Efficiency: $\geq 95\%$
Size: 14.5*9.5cm
Color: White
Instrument Classification: None

01 EPIDEMIC PREVENTION PRODUCTS

3 LAYERS FLAT-FOLD MEDICAL FACE MASK



Elastic earloop with dural point attachment helps provide a secure seal.

Disposable, single-use, breathable, light and environment-friendly.



3 LAYERS CIVIL FACE MASK



Elastic earloop with dural point attachment helps provide a secure seal.

Disposable, single-use, breathable, light and environment-friendly.



PRODUCT PARAMETERS

Specification

Nose Clip: Made of plastic materials
Material: Non-woven Fabric+ Melt-blown Fabric
Sterile: Ethylene oxide disinfection
Flammability: No Afterflame
Supply Ability: 4000000 Piece/Pieces per Day
Standard: GB19083-2010
Certificates: CE

Type: Earloop
Shape: Flat
BEF: ≥95%
Filtration Efficiency: ≥95%
Size: 17.5*9.5cm
Color: Blue/White
Instrument Classification: Class II

PRODUCT PARAMETERS

Specification

Nose Clip: Made of plastic materials
Material: Non-woven Fabric+ Melt-blown Fabric
Sterile: Ethylene oxide disinfection
Flammability: No Afterflame
Supply Ability: 4000000 Piece/Pieces per Day
Standard: GB2626-2006
Certificates: CE

Type: Earloop
Shape: Flat
BEF: ≥95%
Filtration Efficiency: ≥95%
Size: 17.5*9.5cm
Color: Blue/White
Instrument Classification: None

PROTECTIVE COVERALL



ISOLATION GOWN



PRODUCT PARAMETERS

Specification	Packaging Details
Supply Ability: 1000000 Piece/Pieces per Month Leading Time: 20-45 days Shelf Life: 2 years Gender: Unisex Size: S,M,L,XL,XXL,XXXL Standard: EN14126-2003 Place of Origin: China	1 or 5 pc/polybag /PE or opp bag 50 or 100 pcs/ctn
	Certificates
	CE/ISO13485/FDA

PRODUCT PARAMETERS

Specification	Packaging Details
Material: PP,CPE,SMS,SFS,PP+PE Instrument Classification: Class I Color: Blue, White, Green, Yellow etc Supply Ability: 60000 Piece/Pieces per Day Cuff: knitted elastic cuff, rubber elastic cuff Size: S,M,L,XL,XXL,XXXL Place of Origin: China	10PCS/BAG,100PCS/CTN
	Certificates
	CE FDA ISO

NON STERILE SURGEON DOCTOR CAP



ANTI-FOG DISPOSABLE FACE SHIELD



PRODUCT PARAMETERS

Specification

Material: Non-woven Fabric/polypropylene/PP/SMS
Supply Ability: 50000000 Piece/Pieces per Day
Size: 18' 19' 20' 21' 24' or as request
Application: hospital ,medical,clinic
Type: With tie or elastic/ 2ply or 1ply
Instrument Classification: Class I
Certificates: CE , ISO
Feature: Comfortable, disposable,convenient

Shelf Life: 3 years
Gender: Unisex
Color: Blue/Green/White/Yellow
Leading Time: 20-45 days
Weight: 20gsm-40gsm
Packaging Details:
100PC/POLYBAG,
1000PCS/CARTON

PRODUCT PARAMETERS

Specification

Material: PET+Sponge
Supply Ability: 2000000 Piece/Pieces per Day
Size: 33*22cm,32*22cm,32*24cm,28*19cm
Type: Cover all face anti-fog
Instrument Classification: Class I
Certificates: CE, ISO
Feature: Double side anti fog material
Packaging Details: 1pc/bag; 200bag/CTN

Shelf Life: 1 years
Leading Time: 20-45 days
Color: Transparent+blue
Sponge size: 250*30*25mm
Thickness: 0.25mm
Weight: 35g
Standard: GB/T 32610

ANTI-FOG TRANSPARENT PROTECTIVE GLASSES



PRODUCT PARAMETERS

Specification

Material: PVC frame, Anti fog PC lenses
Type: Adjustable Frame, indirect venting
Instrument Classification: Class I
Certificates: CE IS
Feature: Anti UV, Anti fog,
superior scratches and impact resistance
Packaging Details: 1pc/polybag, 10pcs/inner box,
100pcs/ctn. 12 pcs/box, 25boxes/carton; 300 pcs/ctn

Shelf Life: 3 years
Size: Adult
Color: Transparent
Weight: 20 gsm-40 gsm

SHOES COVER BOOT COVER



PRODUCT PARAMETERS

Specification

Type: antistatic shoe cover
Size: 15x36cm, 15x40cm, 17x42cm
Type: Waterproof, With Tie/ or Elastic on the cover
Size: 48*36cm, 38*50cm(Or as per your request)
Supply Ability: 1000000 Piece/Pieces per Day
Weight: 1.4g, 1.8g, 2g, 2.5g, 3g, 4g, 5g, 6g 20-50gsm
Packaging Details: 100pcs/bag, 2000pcs/ctn

Material: Plastic PE, CPE
Color: Blue, White, Green
Shelf Life: 2 years
Sterile: No
Leading Time: 20-45 days
Certificates: CE ISO

LATEX EXAMINATION GLOVE



PRODUCT PARAMETERS

Type	powdered / powder free	Weight	4.5g/pc, 5.0g/pc, 5.5g/pc, 6.0g/pc, 6.5g/pc
Size	XS,S,M,L,XL	Packaging Details	100pcs/box, 10boxes/ctn,3000pcs/bag
Material	nature latex	Certificates	CE ISO

VINYL EXAMINATION GLOVE



PRODUCT PARAMETERS

Type	powdered / powder free	Weight	4.5g/pc, 5.0g/pc, 5.5g/pc
Size	XS,S,M,L,XL	Packaging Details	100pcs/box, 10boxes/ctn,3000pcs/bag
Material	PVC vinyl chloride	Certificates	CE ISO

NITRILE EXAMINATION GLOVE



PRODUCT PARAMETERS

Type	powdered / powder free	Weight	2.5g/pc,3.5g/pc,4.0g/pc, 4.5g/pc,5.0g/pc,5.5g/pc
Size	XS,S,M,L,XL	Packaging Details	100pcs/box, 10boxes/ctn,3000pcs/bag
Material	Nitrile butadiene rubber	Certificates	CE ISO

PE GLOVES



PRODUCT PARAMETERS

Type	powdered / powder free	Weight	As the customers' requirements
Size	Customized	Packaging Details	100pcs/box, 10boxes/ctn,3000pcs/bag
Material	HDPE,LDPE	Certificates	CE ISO