

钢质手动轮椅系列

Fine quality, Integrity and service

Steel Manual Wheelchair Series



Aluminum wheel  
铝合金轮圈

Board table  
餐桌板

Blind rivet  
拉母钉

Flip up armrest  
后掀扶手

Fixed the footrest  
with silk  
螺丝固定脚



KJT601 >>



KJT601-51 >>



KJT603 >>

Fixed armrest and footrest. Chrome-covered frame, PVC cushion. 8" front wheel and 24" rear wheel.  
固定扶手, 固定脚, 车架电镀, PVC 座背垫, 8" 塑胶前轮, 24" 钢丝后轮。

※ 后轮可更换

Fixed armrests, fixed footrest, Chrome-covered frame, PVC back cushion, 8-inch front wheels, 24-inch steel wire rear wheels seat width: 20-inch.  
固定扶手, 固定脚, 车架电镀, PVC 座背垫, 8 寸塑料前轮, 24 寸钢丝后轮, 坐宽 20 寸。

固定扶手, 固定脚, 车架电镀, 牛津布座背垫, 8" 塑胶前轮, 24" 塑料后轮。

※ 后轮可更换

Size(cm)

Model No	座宽	座深	座高	座前缘高	座后缘高	座前缘宽	座后缘宽	座前缘深	座后缘深	座前缘高	座后缘高	重量 (kg)	尺寸 (mm)	重量 (kg)	重量 (kg)	重量 (kg)	重量 (kg)
KJT601	65	24	46	58	20	47	85	102	40	40	100	92.5x22x87	16	18	144	378	
KJT601-51	70	24	51	58	20	47	88	102	40	41	110	93x26x87	18	20	128	316	
KJT603	66	28	46	58	20	47	85	102	40	42	100	92.5x24x88	16	18.5	138	342	



# STEEL MANUAL WHEELCHAIR



**KJT802B >>**

Fixed armrest and footrest. Painted frame. Nylon cushion. 8" Front wheel and 22" rear wheel.  
固定扶手, 固定脚, 车架喷涂, 牛津布坐垫, 8" 塑胶前轮, 22" 钢丝后轮。



# STEEL MANUAL WHEELCHAIR



**KJT605 >>**

Flip up armrest and detachable footrest. Chrome-covered frame, PVC cushion. 8" Front wheel and 24" rear wheel.  
后掀式长扶手, 可拆脚, 车架电镀, PVC座背垫, 8" 塑胶前轮, 24" 钢丝后轮。

※ 后轮可更换



Flip up armrest  
后掀式扶手



Detachable footrest  
可拆脚



**KJT803 >>**

STEEL WHEELCHAIR, detachable backrest & footrest, 8" front castor, 22" rear wheel.  
钢轮椅, 可拆卸靠背和脚板, 8" 前脚轮, 22" 后轮。

※ 后轮可更换



**KJT606 >>**

Detachable armrest and footrest. Chrome-covered frame, PVC cushion. 8" Front wheel and 24" rear wheel.  
可拆短扶手, 可拆脚车架电镀, PVC座背垫, 8" 塑胶前轮, 24" 钢丝后轮。

※ 后轮可更换



Detachable armrest  
可拆扶手



Detachable footrest  
可拆脚

Size(cm)

Model No	Wheelchair	Seat	Backrest	Armrest	Footrest	Front Wheel	Rear Wheel	Seat Height	Seat Depth	Seat Width	Backrest Height	Backrest Depth	Weight (kg)	DIM	NW	GW	CNTS 20' GP FCL	CNTS 40' HQ FCL
KJT802B	62.5	23	46	55	20	48	85	104	40	40	100	96x22x87.5	16	18	146	369		
KJT803	61	20	46	55	20	48	85	104	40	40	100	75x20x79	16	18	228	558		

Size(cm)

Model No	Wheelchair	Seat	Backrest	Armrest	Footrest	Front Wheel	Rear Wheel	Seat Height	Seat Depth	Seat Width	Backrest Height	Backrest Depth	Weight (kg)	DIM	NW	GW	CNTS 20' GP FCL	CNTS 40' HQ FCL
KJT605	63	25	46	58	20	53	91	102	40	42	100	79x28x92	17.5	20	130	324		
KJT606	65	30	46	58	20	47	91	102	40	42	100	79x28x92	17.5	19.5	130	324		





KJT606R >>

Detachable armrests, detachable legs, powder coated frame, 8" front wheels, 24" rear wheels.  
可拆扶手, 可拆脚, 喷涂车架, 8" 前小轮, 24" 后轮。



KJT611 >>

Detachable armrest and footrest, Blue-painted frame. Nylon cushion. 8" PVC front wheel and 24" rear wheel with PU tire.  
可拆手和脚, 蓝色喷涂车架, 牛津布座背垫, 8寸 PVC 前轮, 24寸 PU 后轮, 带背折。



KJT606C >>

Description: Detachable armrests and elevating footrests. Chrome-cover frame, 8-inch front wheels, 24 inch rear wheels.  
可拆扶手, 骨科脚, 电镀车架, 8寸前轮, 24寸后轮。



KJT606-51 >>

Detachable armrests, detachable legs, chromed frame, 8-inch front wheels, 24-inch rear mag wheels.  
可拆扶手, 可拆脚。电镀车架, 8寸前小轮, 24寸塑料后轮



KJT605A >>

Flip up armrest and detachable footrest. Black powdered coated frame, nylon cushion. 8" Front wheel and 24" quick-release wheel.  
后掀式长扶手, 可拆脚, 车架电镀, PVC 座背垫, 8" 塑胶前轮, 24" 快拆后轮。

※ 后轮可更换  
Size(cm)

Model No	Wheelchair	Seat	Backrest	Armrest	Footrest	Seat Height	Seat Width	Seat Depth	Backrest Height	Backrest Width	Weight (Kg)	DIM	NW	GW	CNTS 20" GP FCL	CNTS 40" HQ FCL
KJT606R	64	28	46	58	20	50	91	79	40	40	100	79x28x92	16.7	18.1	130	324
KJT611	65	30	46	60	20	52	85	107	40	41	100	79x28x79	16.5	17.5	154	378
KJT605A	63	25	46	58	20	53	91	102	40	42	100	79x24x92	17.5	20	154	379
KJT605B	63	25	46	58	20	53	91	102	40	42	100	79x28x92	17.5	20	130	324



KJT605B >>

Flip up armrest and detachable footrest. Chrome-covered frame, PVC cushion. 8" Front wheel and 24" rear wheel with IV pole.  
后掀式长扶手, 可拆脚, 车架电镀, PVC 座背垫, 8" 塑胶前轮, 24" 钢丝后轮。

※ 后轮可更换



KJT607 >>

Fixed armrest and footrest with double brakes. Chrome-covered frame, PVC cushion. 8" Front wheel and 24" pneumatic tire rear wheel.  
固定扶手, 固定脚, 车架电镀, PVC 座背垫, 8" 塑胶前轮, 24" 钢丝后轮冲气胎。

※ 后轮可更换  
Size(cm)

Model No	Wheelchair	Seat	Backrest	Armrest	Footrest	Seat Height	Seat Width	Seat Depth	Backrest Height	Backrest Width	Weight (Kg)	DIM	NW	GW	CNTS 20" GP FCL	CNTS 40" HQ FCL
KJT606C	64	28	46	58	20	50	91	79	40	40	100	79x28x92	16.7	18.1	130	324
KJT606-51	71	26	51	58	20	46	86	93	40	44	136	79x28x91	21.5	22.5	119	291
KJT607	65	28	45	60	20	44	86	104	40	42	100	92.5x22x87	16	18	144	378
KJT606B	73	30	51	58	20	50	90	79	40	40	100	81x29x89	18.5	20.5	129	315

KJT606B >>

Flip up armrests, detachable footrest powdered coated frame, 8" front wheels, 24" rear mag wheels.  
后掀式扶手, 可拆脚, 车架喷涂, 8" 前小轮, 24" 塑料后轮





KJT601AM >>

Steel Wheelchair,Frame Color: Powder coating Frame.Fixed Armrest with PVC arm,swing away Footrest,foldable back handle.

铁轮椅,车架颜色:喷粉爆花。固定扶手,可拆脚,折叠背  
※ 后轮可更换



KJT605BM >>

Steel Wheelchair,Frame Color: Powder coating Frame,flip up Armrest with PVC arm,swing away Footrest.

铁轮椅,车架颜色:喷粉爆花,后掀扶手,可拆脚  
※ 后轮可更换



KJT605AM >>

Steel Wheelchair,Frame Color: Powder coating Frame,Flip up Armrest with PVC arm,swing away Footrest.

铁轮椅,车架颜色:蓝色喷涂,后掀扶手,可拆脚  
※ 后轮可更换



KJT605DM >>

Steel Wheelchair,Frame Color: Powder coating Frame,flip up Armrest with PVC arm,rise and fall Footrest,24"rear mag wheel with solid tire

铁轮椅,车架颜色:喷粉爆花,后掀扶手,可拆骨科脚,24"快拆塑料大轮  
※ 后轮可更换

Size(cm)

Model No	Wheelchair	Seat	Back	Wheel	Front	Back	Seat	Back	Seat	Back	Weight	Dimensions	NW	GW	CNTS 20" GP FCL	CNTS 40" HQ FCL
KJT601AM	64	26	46	60	20	52	88	108	40	38	136	78x26x76	17.5	19.5	176	428
KJT605AM	65	28	46	60	20	51	87	108	40	38	136	78x28x88	18	20	140	343

Size(cm)

Model No	Wheelchair	Seat	Back	Wheel	Front	Back	Seat	Back	Seat	Back	Weight	Dimensions	NW	GW	CNTS 20" GP FCL	CNTS 40" HQ FCL
KJT605BM	65	30	46	60	20	52-57	87	109	40	35	136	80x30x88	18.3	20.3	130	315
KJT605DM	65	30	46	60	20	52-57	89	120	40	38	136	80x30x85	20.5	22.5	134	326



# STEEL MANUAL WHEELCHAIR



**KJT802 >>**

Fixed armrest and footrest. Painted frame. Nylon cushion. 8" Front wheel and 24" rear wheel.

固定扶手, 固定脚, 车架喷涂, 牛津格布, 8" 塑胶前轮, 24" 钢丝后轮。

※ 后轮可更换



**KJT804 >>**

Fixed armrest and footrest. Painted frame. PVC cushion. 8" Front wheel and 24" rear wheel.

固定扶手, 固定脚, 车架喷涂, PVC 座背垫, 8" 塑胶前轮, 24" 钢丝后轮。

※ 后轮可更换



**KJT811B >>**

Fixed armrest and footrest. Painted frame. PVC cushion. 8" Front wheel and 24" rear wheel.

固定扶手, 固定脚, 车架喷涂, PVC 座背垫, 8" 塑胶前轮, 24" 钢丝后轮。

※ 后轮可更换



**KJT812 >>**

Fixed armrest and footrest. Painted frame. Nylon cushion. 6" Front wheel and 22" rear wheel.

固定扶手, 固定脚, 车架喷涂, 红格牛津布座背垫, 6" 塑胶前轮, 22" 钢丝后轮。

※ 后轮可更换

Size(cm)

Model No	Wheelchair	Seat	Backrest	Seat Depth	Seat Width	Seat Height	Backrest Height	Backrest Angle	Weight	Weight Capacity	Dimensions	Net Weight	Gross Weight	Front Wheel	Rear Wheel	
KJT802	62.5	25	45	58	20	48	85	104	40	40	100	92.5x22x87	16	18	144	378
KJT804	65	27	45	58	20	43	85	102	40	42	100	92.5x22x87	17	19	144	378
KJT811B	65	27	45	58	20	43	85	102	40	42	100	92.5x22x87	17	19	144	378
KJT812	55	28	38	55	15	42	82	100	40	42	100	89x24x81	14	16	152	381

# STEEL MANUAL WHEELCHAIR



**KJT102 >>**

Fixed armrest and detachable footrest. Painted frame. Nylon cushion. 8" Front wheel and 12" rear wheel.

固定扶手和可拆脚, 喷涂车架, 尼龙座背垫, 8" 前轮, 12" 后轮。



**KJT301-1 >>**

Fixed armrest and detachable footrest. Chrome-cover frame. Nylon cushion. 8" Front wheel and 8" rear wheel.

固定扶手, 可拆脚, 车架电镀, 牛津布座背垫, 8" 前轮, 8" 后轮。

※ 后轮可更换



Size(cm)

Model No	Wheelchair	Seat	Backrest	Seat Depth	Seat Width	Seat Height	Backrest Height	Backrest Angle	Weight	Weight Capacity	Dimensions	Net Weight	Gross Weight	Front Wheel	Rear Wheel	
KJT102	60	26	44	30	20	49	92	94	39	40	136	63x26x79	13	15	208	510
KJT301-1	56	24	46	20	20	46	93	90	40	46	136	63x26x79	11.8	13.8	208	510



Fine quality, Integrity and service

# Commode 座便轮椅系列 wheelchair series

Fashion Design multi-function Commode chair.  
时尚设计, 多功能座便椅。  
The frame can be detachable, so the packing size  
can be small.  
可拆框架, 包装尺寸较小。



KJT728 >>



Back-lifting armrests.  
上掀扶手。



Detachable legs 8" front wheel  
24" quick-release wheel and  
8" rear wheel.  
脚可拆, 8" 前轮, 外加 24" 快拆轮,  
8" 后轮。



Detachable frame.  
可拆框架。

Model No	65	30	44	60	13	54	96	108	41	34	100	85x45x60	16.9	18.9	118	289
KJT728																



KJT704B >>

Detachable long armrests, height adjustable and detachable  
footrest. Retroverted high basket. Seat with commode. Chrome-  
covered frame. 8" PVC front wheel, 24" steel rear wheel.  
可拆长扶手, 可拆式升降脚, 高靠背, 车架电镀, PVC 座便式坐垫, 可卸  
下 .8"PVC 前轮, 24" 钢丝后轮。

※ 后轮可更换

Hydraulic pressure adjust back. Adjust at  
any angle.  
背部采用液压杆调节, 可实现任意角度的调节。



KJT704D >>

Detachable long armrests, detachable lift-able feet, high  
backrest, spraying shelf, PVC cushion for seat, can be  
discharged, 8" front wheel, 24" steel rear wheel.  
可拆长扶手, 可拆式升降脚, 高靠背, 车架电镀, PVC 座便式坐垫, 可  
卸下 .8" 塑胶前轮, 24" 钢丝后轮。

※ 后轮可更换



Model No	63	35	45	58	20	52	120	127	43	78	100	97x32x87	30	34	100	245
KJT704B																
KJT704D																

座便轮椅  
Commode Wheelchair Series



