



Steel frame with powder coating color.
铁制喷粉支架.



Flip up armrest, Easy to transfer users.
上掀式扶手, 方便转移使用者.



Dual-purpose seat plate, free to switch between commode chair or shower chair.
两用座板, 随意切换座便椅或沐浴椅功能.



New design foot pads, to be more stable and slip-resistant.
新型加大脚垫, 使产品更稳固和更防滑.



Soft cushion available.
可选配软坐垫.

KJT739S >>

41-51CM

23CM

Size(cm)

Model No.	Length	Width	Height	N.W.(kg)	G.W.(kg)	D/M	PCS/CN	(kg)	CNTS 20' GP FCL	CNTS 40' HQ FCL
KJT739S	53	53-58	63-73	6	7	55x23x77	1	100	250	650



AS 作为沐浴椅 SHOWER CHAIR

- The shower head can be placed on the clip during bathing.
花洒头通过插座设计可固定在产品旁边, 在沐浴期间使用更方便.
- The cane can be inserted into clip, and can be taken out to assist in standing.
使用者可以把拐杖插入插座中, 在需要站立前可随手拿出, 协助支撑站立.



AS 作为座便椅 COMMODE CHAIR

- Seat cover can be hung on the hook to prevent loss.
产品盖板可悬挂在挂钩上, 防止丢失.
- Put the paper and towel on the hook for easy use.
把纸巾或者手巾挂在挂钩上, 方便使用.



AS 作为马桶座椅 TOILET SEAT

- Hang the seat cover on the hook, and take the bucket away to use as toilet chair.
把产品盖板放在挂钩上, 把便桶取出, 即可作为马桶座椅使用.
- The user can support the product armrest to stand up.
使用者可以支撑产品扶手起立.
- Adjust the height of the product and use it as a toilet seat raiser.
调节产品高度, 可作为马桶增高器使用.



KJT739 >>



KJT740 >>

作为座便椅
COMMODE CHAIR



作为沐浴椅
SHOWER CHAIR



Aluminum frame with matt surface.
哑光铝合金支架。



Stainless steel screws, more rustproof.
不锈钢螺丝，更防锈。



4 level adjustable height.
四级可调高度。

Size(cm)

Model No.	Length	Width	Height	N.W.(kg)	G.W.(kg)	D/M	PCS/CN	(Kg)	CNTS 20' GP FCL	CNTS 40' HQ FCL
KJT739	53	51-56	61-69	3.6	4.6	55x23x77	1	100	290	715
KJT740	53	51-56	70-78	4	5	55x23x88	1	100	240	600



KJT740C >>



Clip design
插座设计

Hook design
挂钩设计



KJT740CS >>



作为沐浴椅
As Shower Chair



作为马桶增高器
As toilet seat raiser

Size(cm)

Model No.	Length	Width	Height	N.W.(kg)	G.W.(kg)	D/M	PCS/CN	(Kg)	CNTS 20' GP FCL	CNTS 40' HQ FCL
KJT740C	53	51-56	70-78	4.6	5.6	55x23x88	1	100	240	600
KJT740CS	53	51-57	70-78	5.8	6.8	56x24x88	1	100	230	560

Commode 座便椅系列 chair Series

The frame can be detachable, and the commode can be fold up. So packing size is small.
框架可拆卸，坐便桶可折叠，因此包装尺寸较小。



Soft Seat Cushion, with detachable cover.
软垫，带可拆盖。



3' Silent Wheels, two with brake, two without brake.
3' 无声车轮，两个带刹，两个不带刹。



Folding footrest, more intimate.
折叠脚踏板，节省空间。



KJT727B >>>

Size(cm)										
Model No.	Length	Width	Height	N.W.(kg)	G.W.(kg)	D/M	PCS/CN	(kg)	CNTS 20' GP FCL	CNTS 40' HQ FCL
KJT727B	65	55	85	6.6	7.6	48x23.5x70	1	100	348	836



Aluminum commode chair with wheels
铁制带轮座便椅
· Light silver aluminum frame/ 铝合金亮银车架
· Padded hard seat/ 加厚硬座板
· Detachable footrest/ 可折叠脚垫

KJT727 >>>

Size(cm)						
Model No	Length	Width	Height	N.W.(kg)	DIM	PCS CN
KJT727	56	45	90	7	53x23.5x64	1



KJT935BH >>>

Multifunctional aluminum shower-commode chair
多功能铝合金沐浴座便椅
· 6 functions in 1/ 一个产品，六个功能
· One button for folding/ 一键收合
· With shower seat and commode/ 带座板和便桶

Size(cm)						
Model No	Length	Width	Height	N.W.(kg)	DIM	PCS CN
KJT935BH	53	45	76-96.3	4.8	54x16x79	1



As Commode chair



As Shower chair

KJT738 >>>

Multifunctional aluminum shower-commode chair
多功能铝合金沐浴座便椅
· 2 functions in 1/ 一个产品，两个功能
· Save space, closing footprint is only 15.5cm/ 节省空间，收合占地面积仅 15.5cm
· Four colors seat board available/ 四色座板可选

Size(cm)						
Model No	Length	Width	Height	N.W.(kg)	DIM	PCS CN
KJT738	48	43	47	2.7	49x29x85	3



KJT-SS715 >>>

**Aluminum folding commode chair
铝合金折叠座便椅**

- Padded PVC soft seat/PVC 软坐垫
- Foldable structure/ 可折叠结构
- Non-slip rubber foot pad/ 防滑脚垫
- With security lock/ 带安全锁
- 3 level adjustable height/ 三级可调高度



KJT715 >>>

**Aluminum folding commode chair
铝合金折叠座便椅**

- Padded PVC hard seat/PVC 硬座板
- Foldable structure/ 可折叠结构
- Non-slip rubber foot pad/ 防滑脚垫
- With security lock/ 带安全锁
- 3 level adjustable height/ 三级可调高度



KJT731 >>>



Heavy-duty Commode Chair, 300 lbs.(135KG)
weight capacity.
加重型座便椅, 可承重 300 磅 (135 公斤)。



KJT-SS716 >>>

**Aluminum folding commode chair
铝合金折叠座便椅**

- Padded PVC soft seat/PVC 软坐垫
- Foldable structure/ 可折叠结构
- Non-slip rubber foot pad/ 防滑脚垫
- 5 level adjustable height/ 五级可调高度



KJT716 >>>

**Aluminum folding commode chair
铝合金折叠座便椅**

- Padded PVC soft seat/PVC 硬坐垫
- Foldable structure/ 可折叠结构
- Non-slip rubber foot pad/ 防滑脚垫
- 5 level adjustable height/ 五级可调高度



KJT732 >>>

Heavy-duty Commode Chair, 650 lbs.(295KG)
weight capacity.
加重型座便椅, 可承重 650 磅 (295 公斤)。



KJT733 >>>

Heavy-duty Commode Chair, 900 lbs.(400KG)
weight capacity.
加重型座便椅, 可承重 900 磅 (400 公斤)。

Size(cm)

Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT-SS715	56	57	81	5.5	58x21x83	1
KJT715	56	57	81	5.5	58x21x83	1

Size(cm)

Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT-SS716	52	55	79.5	4.9	55x18x65	1
KJT716	52	55	79.5	4.9	55x18x65	1

Size(cm)

Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT731	50	60	65	8.2	51x22x65	1
KJT732	52	65	75	7.8	71x54x67	1
KJT733	58	68	75	9.6	78x52x50	1



KJT710B >>



Steel folding commode chair
铁制折叠座便椅



KJT710C >>

Steel folding commode chair
铁制折叠座便椅



KJT710BU >>

Steel folding commode chair
铁制折叠座便椅

Size(cm)						
Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT710B	52	50	75	5	78x41x54	4
KJT710C	52	50	75	5	78x41x54	4

Size(cm)						
Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT710BU	49	52	76	5	78x41x54	4



KJT709 >>

Steel folding commode chair
铁制折叠座便椅



KJT709B >>

Steel folding commode chair
铁制折叠座便椅



KJT709C >>

Steel folding commode chair
铁制折叠座便椅



KJT709D >>

Steel folding commode chair
铁制折叠座便椅



KJT709DU >>

Steel folding commode chair
铁制折叠座便椅



KJT718 >>

Steel folding commode chair
铁制折叠座便椅

Size(cm)						
Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT709	35	36	35	1	70x42x38	10
KJT709B	40	38	38	1.5	73x41x40	8
KJT709C	40	38	46	4	60x51x41	6
KJT709D	42	39	46	2.4	61x57x41	4



KJT713 >>

Steel folding commode chair
铁制折叠座便椅

Size(cm)						
Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT709DU	42	39	46	2.4	61x57x41	4
KJT718	60	50	42	5.73	81x42x59	3
KJT713	43	36	75	2.56	68x58x37	6



KJT705 >>

Steel folding commode chair
铁制折叠座便椅

- Chrome steel frame/ 铁制电镀支架
- 5 level adjustable height/ 五级可调高度
- PE hard seat/PE 硬座板
- With 3" wheels/ 带 3 寸小轮



KJT708 >>

Steel folding commode chair
铁制折叠座便椅

- Chrome steel frame/ 铁制电镀支架
- 5 level adjustable height/ 五级可调高度
- PE hard seat/PE 硬座板



KJT708B >>

Steel folding commode chair
铁制折叠座便椅

- Powder coated steel frame/ 铁制喷涂支架
- 5 level adjustable height/ 五级可调高度
- PE hard seat/PE 硬座板



KJT729 >>

Steel commode chair
铁制座便椅

Size(cm)

Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT729	51	49	71-81	5.5	59x41x42	1



KJT706B >>

Steel commode chair
铁制座便椅

Size(cm)

Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT706B	55	51	77	8	57x28x80	1



KJT720 >>

Aluminum folding commode chair
铝制折叠座便椅

- Matt aluminum frame/ 哑光铝合金支架
- 5 level adjustable height/ 五级可调高度
- PE hard seat/PE 硬座板
- With 3" wheels/ 带 3 寸小轮

Size(cm)

Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT705	53	53	73-84	6.96	50x19x75	1
KJT708	51.5	53	75-85	6.1	50x19x75	1
KJT708B	51.5	53	78-85	6	50x19x75	1



KJT708L >>

Aluminum folding commode chair
铝制折叠座便椅

- Matt aluminum frame/ 哑光铝合金支架
- 5 level adjustable height/ 五级可调高度
- PE hard seat/PE 硬座板

Size(cm)

Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT708L	51	55	47-57	4	51x21x78	1
KJT720	49	55	77	5.2	51x21x78	1



KJT703 >>

Steel commode chair
铁制座便椅

Size(cm)

Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT703	55	53	79	8.2	56x25x64	1



KJT706 >>

Steel commode chair
铁制座便椅

Size(cm)

Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT706	55	52	79	8	51x28x82	1



KJT707C >>>



Steel commode chair with wheels
铁制带轮座便椅

- Powder coated steel frame/ 铁制喷涂车架
- Padded soft seat/ 加厚软坐垫
- Easy to release arm for patient transfers/ 扶手易于释放, 以便转椅病人
- Detachable footrest/ 可拆脚
- With 5" wheels/ 带 5 寸小轮



KJT707 >>>



• can change chrome frame

Steel commode chair with wheels
铁制带轮座便椅

- Blue printed steel frame/ 铁制蓝色喷涂车架
- Padded soft seat/ 加厚软坐垫
- Detachable armrest/ 可拆手
- Detachable footrest/ 可拆脚
- With 5" wheels/ 带五寸小轮

Size(cm)

Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT707C	52	46	93	12	56x54x52	1
KJT707	83	53	85	12	55x55x49	1



KJT-SS717 >>>



Steel commode chair
铁制座便椅

- Powder coated steel frame/ 铁制喷涂支架
- Detachable armrest and footrest/ 可拆手可拆脚
- 7 levels adjustable height/ 七级可调高度



KJT721 >>>



Steel commode chair
铁制座便椅

- Powder coated steel frame/ 铁制喷涂支架
- Fixed armrest and footrest/ 固定手和固定脚
- 7 levels adjustable height/ 七级可调高度

Size(cm)

Model No	Length	Width	Height	N.W. (kg)	DIM	PCS CN
KJT-SS717	53.5	56.5	60-74.5	6	56x24x37	1
KJT721	46	64	62-77	6	63.5x63.5x63.5	4