

Mesh Nebulizer

Item# C101R9



Features

- ▲ Fine particles: large fog, fine atomized particles, particles can reach 1-5 m, drug absorption more fully
- ▲ High speed: microporous spray tablets, atomization speed 0.2ml/min
- ▲ Quieter: with piezoelectric components, quieter design with low noise
- ▲ Flashing light: green light means normal operation, red light means no power or insufficient voltage
- ▲ Energy-saving and environmental protection: skin-friendly and environment-friendly, the mask adopts PP (polypropylene) as the main material, without DEHP (plasticizer/plasticizer)

Speceification	oscillation frequency:100khz±10% condensation rate: ≥0.25ml/min atomized particles: ≤5um medication flow rate:approx.0.4ml/min volume:max.10ml housing material:ABS operating conditions:temperature:+5℃--+40℃ relative humidity: < 80%non-con-densing
Mains connection	10-240V,50-60HZ,0.15A
power	DC5V/1000mA USB power adapter or 2* 1.5v (2 sets of 5 AA alkaline, no battery included)
Dimension & Weight	45(w)*45(h)*100(L)mm 640*410*320mm/50PCS/14KG 278g(without batteries)